SPECIMEN SUBMISSION GUIDELINES

HIV Diagnostics and Reference Laboratory, Diagnostics and Countermeasures Branch Center for Infectious Disease Research Walter Reed Army Institute of Research 508 Research Dr., Silver Spring, MD, 20910

Molecular HIV-1 Viral Load (Quant)/Drug Resistance Test Request Form

Test Requested	Specimen Requirement	Draw Tube	Shipping Conditions
☐ Quant HIV-1 (Viral Load)	☐ 3 ml plasma (PPT or EDTA)	PPT tubes – Invert 8-10X. Spin tubes within 2 hrs. of collection. Centrifuge ≤1300 RCF for 10 minutes. Freeze plasma aliquot at less than -20°C. EDTA tubes – Store blood at 25°C until centrifuged. Centrifuge ≤1300 RCF for 10 minutes within 4 hours of blood collection. Freeze plasma aliquot at less than -20°C.	Frozen -20°C – Ship frozen aliquoted specimen with dry ice if specimen will be received at HDRL after 24 hours of collection. When shipping frozen, use 2+ lbs. dry ice per day of transport. Shipment with an additional 6 lbs. of dry ice is recommended in case of shipment delay.
☐ Sentosa® SQ HIV-1 Genotyping/ Integrase Assay	☐ 3 ml plasma (EDTA only) Viral load MUST BE ≥ 1000 Copies/ml and result must have been obtained within the past 30 days.	EDTA: Invert 8-10 times. Store the tubes upright at room temperature, spin tubes within 2 hours of collection. Centrifuge at 1900 RCF for 10 minutes at 2 to 8°C to remove plasma. Freeze plasma aliquot at -20°C. Viral Load Date Performed	

Please fill the request form completely to ensure timely specimen processing.

PATIENT IDENTIFICATION	CONTACT INFORMATION	CONTACT INFORMATION		
Patient identifiers <u>MUST INCLUDE</u> :	POC			
Full NameOR FMP/SSN	Clinic / Center Center Address	Physician Name Clinic / Center Center Address		
DOB Specimen Draw Date / Time: Ship Date:	Telephone Number Fax Number (Commercial # only; please include area/country cod Alternate POC Name	F. M. william		
PROCESSING LABORATORY (For HDRL use only)				
BARCODE	DATE RECEIVED QUANTITY & TY RECEIVED / INITIA			

Fax/Email a FedEx tracking and/or invoice number to ensure all shipments sent to the HIV Diagnostics and Reference Laboratory are received, IAW CAP GEN.40530
Shipping Address: 508 Research Drive, Silver Spring, MD 20910

Form # TR Aptima Version: March 2025